



# Separation Funnel Shaker by Takahe

With decades of experience in the instrumentation industry, Takahe Analytical Instruments is a leading and trusted name in the field of analytical instrumentation in India. Backed by technological expertise and a rich service experience, our company is a manufacturer of a range of sample preparation instruments from Nitrogen Evaporators, SPEs, Vortexers, and other instruments. Takahe's latest addition to this range is the Separation Funnel Shaker with specifications that meet your requirements.

## Model name: Separation Funnel Shaker

|                              |  |
|------------------------------|--|
| Max no. of samples           | : 6  |
| Sample Volume                | : 50 ml to 250 ml                            |
| Max no. of samples / station | : 3  |
| Stroke per Minute            | : 15 to 150                                  |
| Time                         | : 1 to 999 Minutes                           |
| Pulsed                       | : ON   OFF Times can be set 0 to 999 Seconds |
| Stroke Angle                 | : Adjustable in 5 steps Mechanically         |

## Other details:

|            |                                  |
|------------|----------------------------------|
| Control    | : Microprocessor                 |
| Keys       | : Tactile Switches               |
| BLDC Motor | : Noiseless and maintenance free |
| Wattage    | : 120Watts                       |
| Voltage    | : 110 to 220 Volts               |

\*No of Samples and Station depends on the Sample volume requested and can differ from that mentioned above.

**Contact us at:** B 104, Punit Industrial Premises, Plot No. 11 & 11A,  
TTC Ind. Area, MIDC, Turbhe, Navi Mumbai 400703  
+91 22 2087 4443 | +91 22 2087 9852 | +91 99201 76611  
Email: [info@takaheinstruments.com](mailto:info@takaheinstruments.com) | [www.takaheinstruments.com](http://www.takaheinstruments.com)

